



Institute of Food Research

Protocol from the *Salmonella* group



RNA extraction and purification

- 1) Harvest 2.0 OD₆₀₀ units of bacterial culture (for example 4 mls of a culture with an OD₆₀₀ = 0.5)
- 2) Prepare a 50 ml Falcon tube containing 2/5 of the culture volume ice-cold 5% (v/v) phenol pH 4.3, 95% (v/v) ethanol (Phenol Sigma order code: P-4682). Add the appropriate volume of culture e.g. add 4 mls of culture to 1.6 mls of phenol/ethanol. Stand on ice for at least 30 min but no longer than 2 hours to stabilise the RNA and prevent degradation.
- 3) Centrifuge samples at 3220 X g 4°C for 10 min. Discard supernatant, resuspend bacterial pellets using residual liquid in tubes and transfer to 1.5 ml microcentrifuge tubes.
- 4) Spin tubes 60 secs at maximum speed in a microfuge and discard remaining liquid.
- 5) **Optional:** Freeze pellets at -80°C. Pellets can be kept for up to 1 month before continuing with RNA prep.
- 6) Resuspend pellets in 100 µl TE buffer containing 50 mg/ml lysozyme. Incubate at room temperature for 5 min.
- 7) Add 75 µl lysis reagent (Promega SV Total RNA Purification kit. Cat: Z3100) and mix by inversion several times.
- 8) Add 350 µl RNA dilution buffer from kit (Promega SV, Cat: Z3100). Mix well by inversion.
- 9) Heat samples at 70 °C for 3 min and then centrifuge for 10 min at full speed (13000 rpm).
- 10) Transfer supernatant to clean tubes supplied with the kit. Add 200 µl ethanol and mix by pipette. Transfer to spin columns (Cat: Z3100) and centrifuge columns for 30 secs at full speed. Discard eluate.
- 11) Wash columns with 600 µl wash buffer from the kit. Spin 30 secs at full speed.
- 12) Prepare DNase mix (all reagents supplied with the kit):
 - a) 5 µl 90 mM MnCl₂
 - b) 40 µl DNase core buffer
 - c) 5 µl DNase

- 13) Apply 50 µl of DNase mix to column matrix and incubate at room temperature for 15 min.
- 14) Add 200 µl DNase stop mix (kit) and centrifuge 30 secs at full speed.
- 15) Wash columns with 600 µl wash buffer, by centrifugation for 30 secs at full speed. Discard eluate.
- 16) Wash columns with 250 µl wash buffer, by centrifugation for 30 secs at full speed. Discard eluate.
- 17) Transfer columns to sterile microcentrifuge tubes, add 100 µl RNase-free distilled H₂O and allow to stand for 1 min.
- 18) Centrifuge at 4500 X g for 2 min. Discard column.
- 19) Check the RNA concentration using a spectrophotometer. From exponential phase cultures, you should typically obtain 50-60 µg RNA. Stationary phase and near-stationary phase samples will yield less, but nevertheless plenty to use. Quality of RNA should also be checked by the Aligent Bioanalyzer or a similar method to ensure the RNA is not degraded.

Revised 11.2007